

Centrifuges Clinical

Digital Clinical E-C4.5-8.15CD

- LED Display
- DC motor
- Manual door interlock
- Continuous and timed operation up to 99min
- Includes plastic rotor E-CFR-8.15C and adapters



Digital Pro Clinical E-C4.5-8.15CP

- LCD Display
- Brushless DC motor
- Speed and RCF switch
- Electric door interlock
- Automatic internal diagnostics
- Continuous or timed operation up to 99 min
- Includes plastic rotor E-CFR-8.15C and adapters



Model	Display	Max. Speed (rpm)	Max. RCF (g)	Rotor Plastic	Maximum No. Tubes	Maximum Tube Size (ml)
E-C4.5-8.15CD	LED	4500	2490	Fixed Angle	8	15
E-C4.5-8.15CP	LCD	4500	2490	Fixed Angle	8	15

Accessories: Rotors E-C4.5-8.15CD and E-C4.5-8.15CP

Accessory	Rotor Description	Maximum No. Tubes	Maximum Tube Size (ml)
E-CFR-8.15C	Fixed Angle Plastic Rotor	8	15
E-CFR-8.15C	Fixed Angle Plastic Rotor	12	10 / 7 / 3

Accessories: Tube Adapters E-C4.5-8.15CD and E-C4.5-8.15CP

Accessory	Adapter Description	Rotor	Maximum No. Tubes	Maximum Tube Size (ml)
E-CRA-8.15C-U	Tube Adapter for Collection and Vacutainer Tubes, 18 pcs/pk	E-CFR-8.15C	12	10
E-CPLUG-8.15C	Adapter Plugs for Smaller Tubes, 18 pcs/pk	E-CFA-8.15C-U	12	5



Centrifuges Clinical

Digital Pro Multi-Purpose Clinical E-C6-4.100CP

- LCD Display
- Brushless DC motor
- Speed and RCF switch
- Electric dual door interlock
- Imbalance alarm
- Rotor identification
- Automatic internal diagnostics
- Continuous or timed operation up to 99 min
- Memory capacity 9 programmes
- Acceleration and braking settings
- Includes swing out rotor E-CSR-4.100 and E-CRA-4.100 buckets



Model	Display	Max. Speed (rpm)	Max. RCF (g)	Rotor S/Steel	Maximum No. Tubes	Maximum Tube Size (ml)
E-C6-4.100CP	LCD	6000	4300	Swing Out	4	100

Accessories: Rotors E-C6-4.100CP

Accessory	Rotor Description	Maximum No. Tubes	Maximum Tube Size (ml)
E-CSR-4.100	Swing Out Rotor Body Used with Various Buckets, Max Volume 4 x 100ml	4	100
E-CFR-8.50	Fixed Angle Aluminium Rotor and S/Steel Adapters	8 16	50 15
E-CFR-6.50CS	Fixed Angle Bio-Sealing Aluminium Rotor	6 Culture Tubes	50
E-CFR-30.15	Fixed Angle Aluminium Rotor and S/Steel Adapters	30	15
E-CFR-18.5CS	Fixed Angle Bio-Sealing Aluminium Rotor	18 Culture Tubes	5
E-CFR-12.8PCRS	Fixed Angle Bio-Sealing Aluminium Rotor	12	PCR 8 Strip
E-CSR-2WP	Swing Out Rotor S/Steel Body with Aluminium Buckets	2	96 Well Plate 384 Well Plate

Accessories: Tube Buckets E-C6-4.100CP

Accessory	Bucket Description	Rotor
E-CRA-4.100	4 X 100ml or 4 X 85ml or 4 X 50ml or 8 X 15ml or 20 x 3 – 10ml	E-CSR-4.100
E-CRA-4.100-50	8 X 50ml S/Steel	E-CSR-4.100
E-CRA-4.100-15	16 x 10ml (Blood Collection), 12ml (urine Tube), 15ml S/Steel	E-CSR-4.100
E-CRA-4.100-10	36 x 3 – 10ml (Blood Collection) S/Steel	E-CSR-4.100

