



## Microplate Mixer

### Analogue E-OMM-A

- Analogue speed control
- Orbital shaking movement
- Brushless DC motor
- Continuous operation
- Perfect for 96 well plates
- Wide range of optional accessories available
- Includes single microplate clamp, E-MM-MC



Model	Display	Speed Range (rpm) Single Microplate	Speed Range (rpm) Double Microplate	Orbital Ø (mm)
E-OMM-A	N/A	0 - 1500	0 - 1000	4.5

### Accessories: General E-OMM-A

Accessory	Accessory Description
E-MM-UTP	Universal Top Plate Attachment Ø 100mm for use with E-M-TA1 to E-M-TA5
E-MM-MC	Microplate Clamp
E-MM-DMC	Double Microplate Clamps

### Accessories: Foam Tube Adaptor E-OMM-A

Accessory	Description	No. of Holes	Tube Size Ø (mm)
E-M-TA1-48.6	Foam Tube Adaptor	48	6
E-M-TA2-15.10	Foam Tube Adaptor	15	10
E-M-TA3-12.12	Foam Tube Adaptor	12	12
E-M-TA4-8.16	Foam Tube Adaptor	8	16
E-M-TA5-8.20	Foam Tube Adaptor	8	20



# EINS SCI Product List

- Autoclaves**
  - Laboratory
  - Medical
- Centrifuges**
  - Clinical
  - Haematocrit
  - High Speed
  - Micro
  - Refrigerated
- Dry Bath Incubator**
  - Heating
  - Heating / Cooling
  - Shaking
- Hotplates**
- Hotplate Magnetic Stirrers**
  - Analogue, LED and LCD Displays
  - Multi-Channel
- Liquid Handling**
- Magnetic Stirrers**
  - Analogue and LED Display
  - Cell Culture
  - Multi-Channel
- Microplate Mixers**
- Overhead Stirrers**
- Roller Mixers**
- Rotary Evaporators**
- Rotators**
  - Disk Rotators
  - Rotisserie Rotators
- Shakers**
  - Linear Shakers
  - Orbital Shakers
  - Rocking Shakers
  - 3D Shakers
- Thermal Cyclers**
- Ultrasonic Cleaners**
- Vortex Mixers**

