

EINS SCI

LABORATORY EQUIPMENT CATALOGUE



Overhead Stirrers

Digital E-OHS20-D

- LED display of speed
- Keyless chuck
- Continuous operation
- Brushless DC motor
- Microprocessor controlled gradual acceleration
- Includes E-OHS-US Stand and E-CS5/40-SS Crossed Stirrer



Digital Pro E-OHS20-P, E-OHS40-P

- LCD display to set speed and view speed / torque
- Keyless chuck
- Continuous operation
- Brushless DC motor
- RS232 PC Interface for variable monitoring
- Microprocessor controlled gradual acceleration
- Includes E-OHS-US Stand and E-CS5/40-SS Crossed Stirrer



Model	Display	Speed (rpm)	Max. Stirring Capacity H ₂ O (L)	Power Input (W)
E-OHS20-D	LED	50 - 2200	20	60
E-OHS20-P	LCD	50 - 2200	20	60
E-OHS40-P	LCD	50 - 2200	40	120

Accessories: General E-OHS20-D, E-OHS20-P, E-OHS40-P

Accessory	Accessory Description	Dimensions W x D x H (mm)
E-OHS-US	Universal Stand and Base	200 x 310 x 780

Accessories: Overhead Stirrers Bars E-OHS20-D, E-OHS20-P, E-OHS40-P

Accessory	Material	Stirrer Type	Shaft Length (cm)	Stirrer Size (cm)
E-CS5/40-SS	316 Stainless Steel	Crossed	40	5
E-SBS6/40-SS	316 Stainless Steel	Straight Blade	40	6
E-BS6.8/40-SS	316 Stainless Steel	Blade	40	6.8
E-CFS9/40-SS	316 Stainless Steel	Centrifugal	40	9
E-CS6.5/35-T	PTFE Coated	Crossed	35	6.5
E-SBS7/35-T	PTFE Coated	Straight Blade	35	7
E-BS6.8/35-T	PTFE Coated	Blade	35	6.8
E-CFS8.5/35-T	PTFE Coated	Centrifugal	35	8.5



EINS SCI Product List

- Autoclaves**
 - Laboratory
 - Medical
- Centrifuges**
 - Clinical
 - Haematocrit
 - High Speed
 - Micro
 - Refrigerated
- Dry Bath Incubator**
 - Heating
 - Heating / Cooling
 - Shaking
- Hotplates**
- Hotplate Magnetic Stirrers**
 - Analogue, LED and LCD Displays
 - Multi-Channel
- Liquid Handling**
- Magnetic Stirrers**
 - Analogue and LED Display
 - Cell Culture
 - Multi-Channel
- Microplate Mixers**
- Overhead Stirrers**
- Roller Mixers**
- Rotary Evaporators**
- Rotators**
 - Disk Rotators
 - Rotisserie Rotators
- Shakers**
 - Linear Shakers
 - Orbital Shakers
 - Rocking Shakers
 - 3D Shakers
- Thermal Cyclers**
- Ultrasonic Cleaners**
- Vortex Mixers**

