

EINS SCI

LABORATORY EQUIPMENT CATALOGUE



Orbital Shakers

Analogue E-OS3-A

- Analogue speed control
- Orbital motion
- DC motor
- Continuous operation
- Includes dish attachment and non-slip mat E-S-RMP-1



Digital E-OS3-D

- Double LED independent display for speed and timer
- Digital speed control
- Orbital motion
- DC motor
- Continuous or timed operation up to 1199min
- Includes dish attachment and non-slip mat E-S-RMP-1




Digital Pro E-OS7.5-P

- Double LCD independent display for speed and timer
- Digital speed control
- Orbital motion
- Brushless DC motor
- Continuous or timed operation up to 1199min
- Includes universal attachment and 4 bars E-S-4BH-2



Model	Display	Max. Weight (kg)	Speed (rpm)	Orbital Ø (mm)	Platform Size W x D (mm)
E-OS3-A	N/A	3	0 - 200	20	240 x 215
E-OS3-D	LED	3	40 - 200	20	240 x 215
E-OS7.5-P	LCD	7.5	100 - 500	10	300 x 300



Shaker Accessories

Accessories: Universal Attachment

Accessory	Accessory Description	Shakers (kg)	Platform Size W x D (mm)
E-S-3BH-1	Universal Platform with 3 Horizontal adj. Bars	3	240 x 240
E-S-FCA-1	Universal Platform used with Flask Clamps	3	240 x 240
E-S-3BV-1	Universal Platform with 3 Vertical adj. Bars	3	240 x 240
E-S-RMP-1	Dish Attachment with Non-Slip Mat	3	240 x 215
E-S-BAR-1	Bar for E-S-3BH-1 and E-S-3BV-1	3	N/A
E-S-4BH-2	Universal Platform with 4 Horizontal adj. Bars	5, 7.5 and 10	300 x 300
E-S-FCA-2	Universal Platform used with Flask Clamps	5, 7.5 and 10	320 x 320
E-S-4BV-2	Universal Platform with 4 Vertical adj. Bars	5, 7.5 and 10	300 x 300
E-S-RMP-2	Dish Attachment with Non-Slip Mat	5, 7.5 and 10	300 x 285
E-S-BAR-2	Bar for E-S-4BH-2 and E-S-4BV-2	5, 7.5 and 10	N/A

Accessories: Flask Fixing Clamp Attachment For E-S-FCA-1, E-S-FCA-2

Accessory	Accessory Description	Flask Size (ml)
E-S-FC-25	Flask Fixing Clamp used with E-S-FCA-1 AND E-S-FCA-2	25
E-S-FC-50	Flask Fixing Clamp used with E-S-FCA-1 AND E-S-FCA-2	50
E-S-FC-100	Flask Fixing Clamp used with E-S-FCA-1 AND E-S-FCA-2	100
E-S-FC-200/250	Flask Fixing Clamp used with E-S-FCA-1 AND E-S-FCA-2	200/250
E-S-FC-500	Flask Fixing Clamp used with E-S-FCA-1 AND E-S-FCA-2	500



EINS SCI Product List

- Autoclaves**
 - Laboratory
 - Medical
- Centrifuges**
 - Clinical
 - Haematocrit
 - High Speed
 - Micro
 - Refrigerated
- Dry Bath Incubator**
 - Heating
 - Heating / Cooling
 - Shaking
- Hotplates**
- Hotplate Magnetic Stirrers**
 - Analogue, LED and LCD Displays
 - Multi-Channel
- Liquid Handling**
- Magnetic Stirrers**
 - Analogue and LED Display
 - Cell Culture
 - Multi-Channel
- Microplate Mixers**
- Overhead Stirrers**
- Roller Mixers**
- Rotary Evaporators**
- Rotators**
 - Disk Rotators
 - Rotisserie Rotators
- Shakers**
 - Linear Shakers
 - Orbital Shakers
 - Rocking Shakers
 - 3D Shakers
- Thermal Cyclers**
- Ultrasonic Cleaners**
- Vortex Mixers**

