


**Sliding Lock**  
provides secure storage with easily identified lock and unlock symbols.

**NEW!**

**Slide Rail Design**  
ensures a secure fit allowing easy placement and removal of slides



**Centrally Positioned Numbering**  
enables simple slide organization with easy viewing of location numbers

## 100 Place Premium Plus Slide Box

### Modern design for maximum functionality

- Made from ABS and Polypropylene makes the boxes robust and eliminates potential rusting of metal parts
- Inventory card is removable and replaceable for trouble-free writing on flat surfaces.
- Hinged lid smoothly opens 180° sitting flat on surfaces and prevents misplacement
- EVA foam lines the bottom of the box helping to hold slides securely in place and protect against damage
- Multiple labeling areas, a large area on the lid of the box with "project" and "date" molded in and the back of the index card can be used as the lid is transparent
- Vibrant translucent colors are easy to identify and enable coding options
- Stylish wave surface texture is easy to hold, increasing safety of sample transport
- Stackable, molded feet allow secure horizontal stacking of multiple units for safe transport or storing
- Flexible storage, the units can be stored vertically or horizontally helping to efficiently use space available
- Unique side storage can accommodate a pen, desiccant or other small items needed for working
- Holds 100 - 26 x 76 mm slides
- Patent Pending

### 100-Place Premium Plus Slide Box

Item No.	Color	Item No.	Color
120576	Blue	120579	Red
120577	Orange	120592	Purple
120578	Green	120593	Assorted*

\*Assorted colors (Purple, Blue, Green, Red, Orange)

Lining	Hinge	L x W x H in	L x W x H cm	UOM
EVA Foam	Sliding Latch	8.2 x 7.7 x 1.3	20.7 x 19.5 x 3.2	5 ea/pk 4 pks/case (shrink wrapped packs)

### Inventory Card

Item No.	L x H in	L x H cm	UOM
120618	8.5 x 11	21.6 x 28	10 ea/pk (Polybag with Label) 50 pks/case