



## Vortex Mixers

### Analogue E-VM-A

- Analogue speed control
- Orbital shaking movement
- Shaded pole motor
- Continuous or touch operation
- Constant speed under varying loads
- Silicone feet for ultra-stability
- Wide range of optional accessories available
- Includes standard top cup, E-VM-STC



Model	Display	Speed Range (rpm)	Timer (min)	Orbital Ø (mm)	Outer Dimensions W x H x D (mm)
E-VM-A	N/A	0 - 2500	N/A	4	127 x 130 x 160

### Accessories: General E-VM-A

Accessory	Accessory Description
E-VM-STC	Standard Top Cup for <math>\lt; \text{Ø}30 \text{ mm}</math> tubes
E-VM-UTP	Universal Top Plate with 4 x Rubber Suction Feet (E-VM-RSF), for use with E-M-TA1 to E-M-TA5 or E-VM-PP99
E-VM-PP99	Ø 99 mm Platform Pad for large tubes or small flasks, used with E-VM-UTP
E-VM-RSF	Rubber Suction Foot with Screw, 1pc
E-VM-THR	Tube Holding Rod, for use with E-M-TA1 to E-M-TA5

### Accessories: Foam Tube Adaptor E-VM-A

Accessory	Description	No. of Holes	Tube Size Ø (mm)
E-M-TA1-48.6	Foam Tube Adaptor	48	6
E-M-TA2-15.10	Foam Tube Adaptor	15	10
E-M-TA3-12.12	Foam Tube Adaptor	12	12
E-M-TA4-8.16	Foam Tube Adaptor	8	16
E-M-TA5-8.20	Foam Tube Adaptor	8	20



# EINS SCI Product List

- Autoclaves**
  - Laboratory
  - Medical
- Centrifuges**
  - Clinical
  - Haematocrit
  - High Speed
  - Micro
  - Refrigerated
- Dry Bath Incubator**
  - Heating
  - Heating / Cooling
  - Shaking
- Hotplates**
- Hotplate Magnetic Stirrers**
  - Analogue, LED and LCD Displays
  - Multi-Channel
- Liquid Handling**
- Magnetic Stirrers**
  - Analogue and LED Display
  - Cell Culture
  - Multi-Channel
- Microplate Mixers**
- Overhead Stirrers**
- Roller Mixers**
- Rotary Evaporators**
- Rotators**
  - Disk Rotators
  - Rotisserie Rotators
- Shakers**
  - Linear Shakers
  - Orbital Shakers
  - Rocking Shakers
  - 3D Shakers
- Thermal Cyclers**
- Ultrasonic Cleaners**
- Vortex Mixers**

