

Cardboard Cryogenic Storage Boxes

- Store in mechanical freezers or vapor phase LN₂, -196°C. Must use tubes capable of being stored at the required temperature of the user.
- Not recommended for autoclaving.

Cryogenic Vial Boxes & Partitions

For storage of small samples and tubes

- Made from cardboard
- Colors allow for easy identification and sorting purposes
- All box footprints are 13.3 x 13.3 cm, select from two different heights



Cardboard Cryogenic Vial Boxes

Partitions must be ordered separately

Item No.	Color	Tubes mL	L x W x H in	L x W x H cm	UOM
HS2860A	White	<input type="checkbox"/> 1.5/2.0	5.2 x 5.2 x 2	13.3 x 13.3 x 5	12 ea/pk, 4 pk/cs
HS2860B	White	<input type="checkbox"/>			
HS2860C	Assorted	<input type="checkbox"/>			
HS2860CB	Blue	<input checked="" type="checkbox"/> 1.5/2.0	5.2 x 5.2 x 3	13.3 x 13.3 x 7.5	12 ea/pk, 4 pk/cs
HS2860CR	Red	<input checked="" type="checkbox"/>			
HS2860CY	Yellow	<input checked="" type="checkbox"/>			



Vent Holes

allow faster cooling and removal of LN₂

Partitions

Item No.	Array	Tubes mL	Wells	Well Dia. mm	UOM
HS2860E	10 x 10	1.5/2.0	100	11.5	
HS2860H	9 x 9	1.5/2.0	81	12	12 ea/pk
HS2860F	8 x 8	1.5/2.0	64	15	40 pk/cs
HS2860G	7 x 7	1.5/2.0	49	17.5	

Microtube Storage Boxes

Safe microtube storage

Low profile boxes are available in two sizes to accommodate standard 0.5 mL and 0.2 mL microtubes and vessels of a similar size.


- Partitions included in the boxes
- Will fit into most standard freezer racking



Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	UOM
HS120220	White	<input type="checkbox"/> 0.2	169	9.5	5.2 x 5.2 x 1.1	13.3 x 13.3 x 2.9	10 ea/pk
HS120221	White	<input type="checkbox"/> 0.5	100	12	5.2 x 5.2 x 1.6	13.3 x 13.3 x 4.1	4 pk/cs

15/50 mL Tube Storage Boxes & Partitions

- Economical cardboard boxes are disposable
- Write-on surface or descriptions for self numbering systems
- Partitions available in two sizes



Cardboard Tube Storage Box

Partitions must be ordered separately

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	UOM
HS2860M	White	<input type="checkbox"/>					
Box & Lid	15/50	—	5.7 x 5.7 x 4.8	14.5 x 14.5 x 12.2	10 ea/pk, 4 pk/cs		

Partitions

Item No.	Array	Tubes mL	Wells	Well Dia. mm	UOM
HS2860T	6 x 6	15	36	17.5	10 ea/pk, 20 pk/cs
HS2860P	4 x 4	50	16	29.5	

Cryogenic Vial Tube Storage Box

- Two-piece set-up boxes with finger holds
- Boxes made from paper coated cardboard
- 5x5 format for tubes up to 14mm in diameter
- Partitions are included in each box



NEW!

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H mm	UOM
120519	Blue	<input checked="" type="checkbox"/> 2.0	25	14	3.4 x 3.4 x 1.5	87x87x38	25 ea/pk 8 pk/cs

Related Products:

Cryogenic Vial Rack	30
True North® Cryogenic Vials	47
Micro & Conical Centrifuge Tubes	47
Freezer Racks	50