

# EINS SCI


## LABORATORY EQUIPMENT CATALOGUE



## Rotisserie Rotators

### Analogue E-RRM-A

- Analogue speed control
- Continuous operation
- DC motor
- Suitable for cold rooms and incubators between 5 and 40 °C
- Includes 36 x 1.5ml, 12 x 15ml and 12 x 50ml plastic tube clips



Model	Display	Speed Range (rpm)	Timer (min)	Outer Dimensions W x H x D (mm)
E-RRM-A	N/A	0 - 80	N/A	530 x 190 x 150

### Accessories: E-RRM-A

Accessory	Accessory Description	No. of Clips	Tube Size (ml)
E-RAA-36.1.5	Centrifuge Plastic Clips	36	1.5
E-RAA-12.15	Centrifuge Plastic Clips	12	15
E-RAA-12.50	Centrifuge Plastic Clips	12	50

### Digital Pro E-RRM-P

- LCD display for speed and timer
- Digital speed control
- Continuous or timed operation
- DC motor
- Suitable for cold rooms and incubators between 5 and 40 °C
- Standard with stainless steel rotor and 16 x 50 ml stainless steel clips



Model	Display	Speed Range (rpm)	Timer (min)	Outer Dimensions W x H x D (mm)
E-RRM-P	LCD	10 - 70	1 - 1199	510 x 260 x 220

### Accessories: E-RRM-P

Accessory	Accessory Description	No. of Clips	Tube Size (ml)
E-RRA-48/1.5H	S/Steel Rotor and Clips Horizontal	48	1.5
E-RRA-32/1.5V	S/Steel Rotor and Clips Vertical	32	1.5
E-RRA-24/15H	S/Steel Rotor and Clips Horizontal	24	15
E-RRA-16/15V	S/Steel Rotor and Clips Vertical	16	15
E-RRA-24/50H	S/Steel Rotor and Clips Horizontal	24	50
E-RRA-16/50V	S/Steel Rotor and Clips Vertical	16	50



# EINS SCI Product List

- Autoclaves**
  - Laboratory
  - Medical
- Centrifuges**
  - Clinical
  - Haematocrit
  - High Speed
  - Micro
  - Refrigerated
- Dry Bath Incubator**
  - Heating
  - Heating / Cooling
  - Shaking
- Hotplates**
- Hotplate Magnetic Stirrers**
  - Analogue, LED and LCD Displays
  - Multi-Channel
- Liquid Handling**
- Magnetic Stirrers**
  - Analogue and LED Display
  - Cell Culture
  - Multi-Channel
- Microplate Mixers**
- Overhead Stirrers**
- Roller Mixers**
- Rotary Evaporators**
- Rotators**
  - Disk Rotators
  - Rotisserie Rotators
- Shakers**
  - Linear Shakers
  - Orbital Shakers
  - Rocking Shakers
  - 3D Shakers
- Thermal Cyclers**
- Ultrasonic Cleaners**
- Vortex Mixers**

