

Magnetic Stirrers

Analogue E-MS3-A

- Analogue speed control
- ABS casing
- Raised control panel to avoid liquid spillage on controls
- Maximum Magnetic Bar length 50mm
- Round ABS work-plate has good chemical resistance



Digital E-MS3-D

- LED display for speed reading and control
- ABS casing
- Raised control panel to avoid liquid spillage on controls
- Maximum Magnetic Bar length 50mm
- Round ABS work-plate has good chemical resistance



Model	Display	Max. Stirring (L) H ₂ O	Stirring Speed (rpm)	Work-Plate Ø (mm)
E-MS3-A	N/A	3	0 - 1500	135
E-MS3-D	LED	3	100 - 1500	135

Ultra-Thin E-MS0.8-U

- Push button speed control
- ABS casing
- Automatic reverse rotation for thorough mixing
- Maximum Magnetic Bar length 40mm
- Flat PET work-plate has good chemical resistance



Model	Display	Max. Stirring (L) H ₂ O	Stirring Speed (rpm)	Work-Plate Ø (mm)
E-MS0.8-U	N/A	0.8	15 - 1500	90



Note: All Magnetic Stirrers do NOT come with magnetic stirring bars. Wide range of different sizes available on request.

Magnetic Stirrers

Analogue E-MS20-A

- Analogue speed control
- Cast aluminium casing
- Raised control panel to avoid liquid spillage on controls
- Maximum magnetic bar length 80mm
- Round stainless steel ceramic coated work-plate is extremely chemical resistant



Model	Display	Max. Stirring (L) H ₂ O	Stirring Speed (rpm)	Work-Plate Ø (mm)
E-MS20-A	N/A	20	0 - 1500	135

Accessories: E-MS20-A

Accessory	Accessory Description
E-MSH-PC1	Protective Cover, Temperature Resistant to 135°C

Analogue E-MS10S-A

- Analogue speed control
- Cast aluminium casing
- Raised control panel to avoid liquid spillage on controls
- Maximum magnetic bar length 80mm
- Square solid glass ceramic work-plate is extremely chemical resistant



Model	Display	Max. Stirring (L) H ₂ O	Stirring Speed (rpm)	Work-Plate W x D (mm)
E-MS10S-A	N/A	10	0 - 1500	184 x 184



Note: All Magnetic Stirrers do NOT come with magnetic stirring bars. Wide range of different sizes available on request.

EINS SCI Product List

- Autoclaves**
 - Laboratory
 - Medical
- Centrifuges**
 - Clinical
 - Haematocrit
 - High Speed
 - Micro
 - Refrigerated
- Dry Bath Incubator**
 - Heating
 - Heating / Cooling
 - Shaking
- Hotplates**
- Hotplate Magnetic Stirrers**
 - Analogue, LED and LCD Displays
 - Multi-Channel
- Liquid Handling**
- Magnetic Stirrers**
 - Analogue and LED Display
 - Cell Culture
 - Multi-Channel
- Microplate Mixers**
- Overhead Stirrers**
- Roller Mixers**
- Rotary Evaporators**
- Rotators**
 - Disk Rotators
 - Rotisserie Rotators
- Shakers**
 - Linear Shakers
 - Orbital Shakers
 - Rocking Shakers
 - 3D Shakers
- Thermal Cyclers**
- Ultrasonic Cleaners**
- Vortex Mixers**

